



## CASE REPORT

Mötus AX & Mötus AY

Alexandrite

## HYPERPIGMENTATION

**Patient data:**  
27-year old female patient presenting sun spots on her back.

### Alexandrite Laser Settings

Pulse Duration	15 ms
Fluence	24 J/cm <sup>2</sup>
Frequency	2 Hz
Handpiece	10 mm
Number of Sessions	1



*Before treatment.*



*21 days after 1 single session.*

**DEKA**  
Innate Ability

**DEKA M.E.L.A. s.r.l.**  
Via Baldanzese, 17  
50041 Calenzano (FI) - Italy  
Ph. +39 055 88.74.942  
Fax +39 055 88.32.884



[www.dekalaser.com](http://www.dekalaser.com)